

CE  
2460

# EYECRYL™ PLUS CLEAR

Monofocal Hydrophilic Acrylic Aspheric  
Foldable Intraocular Lens



biotech

VISION CARE

# EYECRYL™ PLUS CLEAR

Monofoveal Hydrophilic Acrylic Aspheric  
Foldable Intraocular Lens


## FEATURES

- Aspheric optic surface
- 360 ° square edge
- Increased hydrophobicity because of surface modification hydrophilic IOL
- Less attractive to cell attachments & cell growth
- Surface modification helps in preventing the secondary cataract
- Exact refractive value of the IOL is printed along with diopter on the box
- 5° angulation ensures better contact with posterior capsule
- UV blocking Material
- Exact refractive value of the IOL is printed along with diopter on the box

## SPECIFICATIONS

MATERIAL	Hydrophilic Acrylic CO UV									
OPTIC TYPE	Single Piece, 360° Square Edge with Aspheric Optic									
OPTIC SIZE	6.00 mm									
OVERALL SIZE	12.50 mm									
ANGULATION	5°									
ACD	4.96									
REFRACTIVE INDEX	1.462									
RECOMMENDED ULTRASOUND A-CONSTANT	SRK-T 118.0									
RECOMMENDED OPTICAL A-CONSTANTS	<table><tr><td>SRK - T 118.3</td><td>SRK - II 118.5</td><td>Holl 1 Const SF : 1.37</td></tr><tr><td colspan="2">HOFFER Q ACD : 5.13</td><td>HAIGIS a0:1.030, a1:0.40, a2:0.10</td></tr><tr><td>Barrett 1.52</td><td>Hill-RBF 118.3</td><td>Olsen ACD 4.16</td></tr></table>	SRK - T 118.3	SRK - II 118.5	Holl 1 Const SF : 1.37	HOFFER Q ACD : 5.13		HAIGIS a0:1.030, a1:0.40, a2:0.10	Barrett 1.52	Hill-RBF 118.3	Olsen ACD 4.16
SRK - T 118.3	SRK - II 118.5	Holl 1 Const SF : 1.37								
HOFFER Q ACD : 5.13		HAIGIS a0:1.030, a1:0.40, a2:0.10								
Barrett 1.52	Hill-RBF 118.3	Olsen ACD 4.16								
DIOPTR RANGE	+5.0 D to +30.0 D ( with 0.5 D steps )									
IMPLANTATION SITE	Capsular Bag									
STERILIZATION	Steam									
SHELF LIFE	4 years from date of manufacture									

## Model: HSAS 600



Scan For Biotech Connect App